




The "X" of the SonoAce X8 represents 'eXtreme quality.'  
The SonoAce X8, the mainstay of the X series, demonstrates  
our will and ability to provide you with superb image quality  
and maximum performance.



Aortic valve regurgitation in color Doppler

### Quick Diagnoses and Quick Operation to Accelerate Patient Throughput

With its large LCD monitor and high frame rate, the SonoAce X8 provides higher image quality with more accuracy and speed. Through various high-end ultrasound technologies, the SonoAce X8 helps you to diagnose quickly and easily.


- High resolution 17" LCD monitor
- Speckle Reduction Filter™
- Full Spectrum Imaging™
- DynamicMR™
- High-sensitivity Powerful Doppler performance
- Dual live mode
- Versatile back-up support - built in DVD RW / DICOM SR / 6 USB Ports / USB 2.0 / External HDD / Fast back-up and restoration export system
- Optimized parameters for 2D and Doppler
- Optimized parameters for real-time 3D and 4D



# Speedy

*The SonoAce X8 provides  
accurate and speedy diagnosis in  
your daily practice.*





Provides a monitor arm with free up and down and right and left movement.

**Articulated monitor arm**



Front and rear handles enable you easy and fast movement.



Ergonomic and innovative design allows optimum use of space.

# Slim

Ergonomic and innovative design allows free movement and optimum use of space.



**17 inches flat LCD monitor** Supports high resolution upto 1280x1024 pixels.



**Rotating control panel**

Enables you to adopt the most convenient posture.

**Front and rear handles**

Front and rear handles provide you with ease of movement.




**Internal cables**


The complicated wiring on the rear has been placed inside the unit.

**4-wheel swivel**

Enables smooth movement.







Umbilical artery in Doppler

**Post Measurement**

Enables you to measure and analyze volume data after completion of scanning.

**Post Processing**

Enhances image quality by optimizing parameters after scanning.

**Pointer Function**

Allows for reviewing and measuring saved images while scanning.

**Quadrant Image**

Displays four images simultaneously on one screen and helps you analyze clinical information in no time.

**QuickScan™**

Provides optimal image brightness and contrast automatically with one click.



**Smart**

*A software-based system structure offers ease of upgrade and makes various high-tech functions possible.*

# Multiple Functionality in all Application Fields

The SonoAce X8 stands out in any environment and delivers powerful performance, convenient operation and slim mobility in any application field.

Speedy
Slim
Smart



#### *Applications covered*

- OB/GYN
- 3D/4D
- General Imaging
- Cardiology